

# UGT1A1 Polyclonal Antibody



Catalog Number: 164759

---

## Product name

UGT1A1 Polyclonal Antibody

## Specificity

Human, Mouse, Rat

## Antibody description

Polyclonal antibody to UGT1A1

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-200 of human UGT1A1 (NP\_000454.1).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


WB, IF

## Dilutions

WB 1:500 - 1:2000

IF 1:50 - 1:200

## Validations



Immunofluorescence - UGT1A1 Polyclonal Antibody

Immunofluorescence analysis of U2OS cells using UGT1A1 antibody . Blue: DAPI for nuclear staining.