

SMYD5 Polyclonal Antibody



Catalog Number: 164763

Product name

SMYD5 Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to SMYD5

Preparation

Antigen: Recombinant protein of human SMYD5

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC, IF

Dilutions

WB 1:500 - 1:2000

IHC 1:50 - 1:100

IF 1:50 - 1:100

Validations



Western blot - SMYD5 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using SMYD5 antibody .Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution.Lysates/proteins: 25ug per lane.Blocking buffer: 3% nonfat dry milk in TBST.

Immunohistochemistry - SMYD5 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human kidney using SMYD5 antibody .

Immunofluorescence - SMYD5 Polyclonal Antibody

Immunofluorescence analysis of A549 cells using SMYD5 antibody . Blue: DAPI for nuclear staining.