

BNIP3L Polyclonal Antibody



Catalog Number: 164798

Product name

BNIP3L Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to BNIP3L

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-187 of human BNIP3L (NP_004322.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IHC

Dilutions

WB 1:500 - 1:2000


IHC 1:50 - 1:200

Validations




Western blot - BNIP3L Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using BNIP3L antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 30s.



Immunohistochemistry - BNIP3L Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human stomach using BNIP3L antibody at dilution



of 1:100 (40x lens).


BNIP3L Polyclonal Antibody



Catalog Number: 164798

Immunohistochemistry - BNIP3L Polyclonal Antibody

Immunohistochemistry of paraffin-embedded rat brain using BNIP3L antibody at dilution of 1:100



(40x lens).

Immunohistochemistry - BNIP3L Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human lung cancer using BNIP3L antibody at dilution of 1:100 (40x lens).