

UBE2V1 Polyclonal Antibody



Catalog Number: 164812

Product name

UBE2V1 Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to UBE2V1

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-147 of human UBE2V1 (NP_001027459.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IF

Dilutions

WB 1:500 - 1:2000


IF 1:10 - 1:100

Validations




Western blot - UBE2V1 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using UBE2V1 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 90s.



Immunofluorescence - UBE2V1 Polyclonal Antibody

Immunofluorescence analysis of U2OS cells using UBE2V1 antibody.



UBE2V1 Polyclonal Antibody



Catalog Number: 164812

Immunofluorescence - UBE2V1 Polyclonal
Antibody

Immunofluorescence analysis of A549 cells using
UBE2V1 antibody .