

# SOAT1 Polyclonal Antibody



Catalog Number: 164813

---

**Product name**

Polyclonal

SOAT1 Polyclonal Antibody

**Ig Type****Specificity**

Rabbit IgG

Human, Mouse, Rat

**Applications****Antibody description**

WB

Polyclonal antibody to SOAT1

**Dilutions****Preparation**

WB 1:500 - 1:2000

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-130 of human SOAT1 (NP\_003092.4).

**Validations**

Western blot - SOAT1 Polyclonal Antibody

**Formulation**

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Western blot analysis of extracts of various cell lines, using SOAT1 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 10s.

**Storage**

Store at -20°C. Avoid freeze / thaw cycles.

**Clonality**