

# C3AR1 Polyclonal Antibody



Catalog Number: 164835

## Product name

C3AR1 Polyclonal Antibody

## Specificity

Human, Mouse

## Antibody description

Polyclonal antibody to C3AR1

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 161-337 of human C3AR1 (NP\_004045.1).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

WB

## Dilutions

WB 1:500 - 1:2000

## Validations



## Western blot - C3AR1 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using C3AR1 Antibody at 1:3000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 15s.