

# FGF9 Polyclonal Antibody



Catalog Number: 164842

## General Information

### Product name:

FGF9 Polyclonal Antibody

### Immunogen:

### Applications:

WB, IHC, IF

### Specificity:

Human, Mouse, Rat

### Antibody description:

Polyclonal antibody to FGF9

## Properties

### Preparation:

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-208 of human FGF9 (NP\_002001.1).

### Formulation:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

### Storage:

Store at -20°C. Avoid freeze / thaw cycles.

### Clonality:

Polyclonal

### Ig Type:Rabbit IgG


## Applications

WB 1:500 - 1:2000

IHC 1:50 - 1:100

IF 1:50 - 1:100

## Validations



### Western blot - FGF9 Polyclonal Antibody

Western blot analysis of extracts of Mouse kidney, using FGF9 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 90s.