

# IREB2 Polyclonal Antibody



Catalog Number: 164845

---

## Product name

IREB2 Polyclonal Antibody

## Specificity

Human

## Antibody description

Polyclonal antibody to IREB2

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-250 of human IREB2 (NP\_004127.1).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

IF

## Dilutions

IF 1:20 - 1:50

## Validations



Immunofluorescence - IREB2 Polyclonal Antibody

Immunofluorescence analysis of A549 cells using IREB2 antibody . Blue: DAPI for nuclear staining.