

TRIM27 Polyclonal Antibody



Catalog Number: 164853

Product name

TRIM27 Polyclonal Antibody

Specificity

Human, Mouse

Antibody description

Polyclonal antibody to TRIM27

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 214-513 of human TRIM27 (NP_006501.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC, IF


Dilutions

WB 1:500 - 1:2000

IHC 1:50 - 1:200


IF 1:10 - 1:100

Validations



Western blot - TRIM27 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using TRIM27 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking



buffer: 3% nonfat dry milk in TBST.


Western blot - TRIM27 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using TRIM27 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST.

TRIM27 Polyclonal Antibody




Catalog Number: 164853




Immunohistochemistry - TRIM27 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human esophagus using TRIM27 Antibody at dilution of 1:100 (40x lens).



Immunohistochemistry - TRIM27 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded mouse stomach using TRIM27 Antibody at dilution



of 1:100 (40x lens).

Immunofluorescence - TRIM27 Polyclonal Antibody

Immunofluorescence analysis of U2OS cells using TRIM27 antibody . Blue: DAPI for nuclear staining.