

# KLF10 Polyclonal Antibody



Catalog Number: 164861

---

## Product name

KLF10 Polyclonal Antibody

## Specificity

Human

## Antibody description

Polyclonal antibody to KLF10

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-250 of human KLF10 (NP\_001027453.1).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


WB, IF

## Dilutions

WB 1:500 - 1:2000

IF 1:50 - 1:100

## Validations



Immunofluorescence - KLF10 Polyclonal Antibody

Immunofluorescence analysis of MCF7 cells using KLF10 antibody . Blue: DAPI for nuclear staining.