

# PRAME Polyclonal Antibody



Catalog Number: 164886

---

## Product name

PRAME Polyclonal Antibody

## Specificity

Human

## Antibody description

Polyclonal antibody to PRAME

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human PRAME (NP\_996839.1).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


WB, IF

## Dilutions

WB 1:500 - 1:2000

IF 1:20 - 1:50

## Validations



Immunofluorescence - PRAME Polyclonal Antibody

Immunofluorescence analysis of U2OS cells using PRAME antibody . Blue: DAPI for nuclear staining.