

RHOD Polyclonal Antibody



Catalog Number: 164889

Product name

RHOD Polyclonal Antibody

Specificity

Human

Antibody description

Polyclonal antibody to RHOD

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-210 of human RHOD (NP_055393.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IF

Dilutions

WB 1:500 - 1:2000

IF 1:50 - 1:100

Validations



Immunofluorescence - RHOD Polyclonal Antibody

Immunofluorescence analysis of MCF-7 cells using RHOD antibody . Blue: DAPI for nuclear staining.