

AHSP Polyclonal Antibody



Catalog Number: 164891

Product name

Polyclonal

AHSP Polyclonal Antibody

Ig Type

Specificity

Rabbit IgG

Human, Mouse, Rat

Applications

Antibody description

WB, IHC

Polyclonal antibody to AHSP

Dilutions

Preparation

WB 1:500 - 1:2000

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-102 of human AHSP (NP_057717.1).

IHC 1:50 - 1:200


Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Validations

Storage

Store at -20°C. Avoid freeze / thaw cycles.



Clonality

Western blot - AHSP Polyclonal Antibody