

TRPV5 Polyclonal Antibody



Catalog Number: 164894

Product name

TRPV5 Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to TRPV5

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-325 of human TRPV5 (NP_062815.3).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IF

Dilutions

WB 1:500 - 1:2000


IF 1:10 - 1:100

Validations



Western blot - TRPV5 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using TRPV5 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 90s.



Immunofluorescence - TRPV5 Polyclonal Antibody

Immunofluorescence analysis of U2OS cells using TRPV5 antibody.