

# POMGNT2 Polyclonal Antibody



Catalog Number: 164897

---

## Product name

POMGNT2 Polyclonal Antibody

## Specificity

Human

## Antibody description

Polyclonal antibody to POMGNT2

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 301-580 of human POMGNT2 (NP\_116195.2).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

IF

## Dilutions

IF 1:20 - 1:50

## Validations



Immunofluorescence - POMGNT2 Polyclonal Antibody

Immunofluorescence analysis of A549 cells using POMGNT2 antibody . Blue: DAPI for nuclear staining.