

LMNB2 Polyclonal Antibody



Catalog Number: 164898

Product name

LMNB2 Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to LMNB2

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 341-620 of human LMNB2 (NP_116126.3).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC, IF


Dilutions

WB 1:500 - 1:2000

IHC 1:50 - 1:200


IF 1:50 - 1:100

Validations




Western blot - LMNB2 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using LMNB2 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 90s.



Immunohistochemistry - LMNB2 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded rat brain using LMNB2 Antibody at dilution of 1:200





(40x lens).

LMNB2 Polyclonal Antibody

Catalog Number: 164898

Immunohistochemistry - LMNB2 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded rat pancreas using LMNB2 Antibody at dilution of 1:200 (40x lens).





Immunohistochemistry - LMNB2 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human gastric cancer using LMNB2 Antibody at dilution of 1:200 (40x lens).

Immunohistochemistry - LMNB2 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human liver injury using LMNB2 Antibody at dilution of 1:200 (40x lens).




Immunohistochemistry - LMNB2 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded mouse brain using LMNB2 Antibody at dilution of

Immunohistochemistry - LMNB2 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human colon carcinoma using LMNB2 Antibody at dilution of 1:200 (40x lens).



1:200 (40x lens).

Immunohistochemistry - LMNB2 Polyclonal Antibody

LMNB2 Polyclonal Antibody



Catalog Number: 164898

Immunohistochemistry of paraffin-embedded rat brain using LMNB2 Antibody at dilution of 1:100 (40x lens).