

# ASMT Polyclonal Antibody



Catalog Number: 164922

---

## Product name

ASMT Polyclonal Antibody

## Specificity

Human

## Antibody description

Polyclonal antibody to ASMT

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-298 of human ASMT (NP\_001164510.1).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


WB, IF

## Dilutions

WB 1:500 - 1:2000

IF 1:50 - 1:100

## Validations



Immunofluorescence - ASMT Polyclonal Antibody

Immunofluorescence analysis of MCF7 cells using ASMT antibody .