

# BAP1 Polyclonal Antibody



Catalog Number: 164923

## Product name

BAP1 Polyclonal Antibody

WB 1:500 - 1:2000

IHC 1:50 - 1:200

IF 1:50 - 1:200

## Specificity

Human, Mouse, Rat

## Antibody description

Polyclonal antibody to BAP1

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 460-729 of human BAP1 (NP\_004647.1).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type


Rabbit IgG

## Applications

WB, IHC, IF


## Dilutions

## Validations



Immunohistochemistry - BAP1 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human breast using BAP1 antibody at dilution of



1:150 (40x lens).

Immunohistochemistry - BAP1 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human breast cancer using BAP1 antibody at dilution of 1:150 (40x lens).