

# HYAL1 Polyclonal Antibody



Catalog Number: 164957

## Product name

HYAL1 Polyclonal Antibody

## Specificity

Human, Mouse, Rat

## Antibody description

Polyclonal antibody to HYAL1

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 136-435 of human HYAL1 (NP\_149349.2).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


WB, IHC

## Dilutions

WB 1:500 - 1:2000


IHC 1:50 - 1:200

## Validations



Western blot - HYAL1 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using HYAL1 Antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 90s.



Immunohistochemistry - HYAL1 Polyclonal Antibody


Immunohistochemistry of paraffin-embedded human liver injury using HYAL1 antibody at dilution of 1:100 (40x lens).

# HYAL1 Polyclonal Antibody




Catalog Number: 164957

---



Immunohistochemistry - HYAL1 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human stomach using HYAL1 antibody at dilution



of 1:100 (40x lens).

Immunohistochemistry - HYAL1 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded rat kidney using HYAL1 antibody at dilution of 1:100 (40x lens).