

TAF1C Polyclonal Antibody



Catalog Number: 165026

Product name

TAF1C Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to TAF1C

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 70-320 of human TAF1C (NP_001230085.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IHC

Dilutions

WB 1:500 - 1:2000


IHC 1:50 - 1:200

Validations




Western blot - TAF1C Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using TAF1C antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 90s.



Immunohistochemistry - TAF1C Polyclonal Antibody

Immunohistochemistry of paraffin-embedded mouse lung using TAF1C Antibody at dilution of



1:100 (40x lens).


TAF1C Polyclonal Antibody



Catalog Number: 165026

Immunohistochemistry - TAF1C Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human gastric cancer using TAF1C Antibody at dilution of 1:100 (40x lens).



Immunohistochemistry - TAF1C Polyclonal Antibody

Immunohistochemistry of paraffin-embedded mouse lung using TAF1C Antibody at dilution of 1:100 (40x lens).