

# TEAD1 Polyclonal Antibody



Catalog Number: 165027

## Product name

TEAD1 Polyclonal Antibody

## Specificity

Human, Mouse, Rat

## Antibody description

Polyclonal antibody to TEAD1

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 135-215 of human TEAD1 (NP\_068780.2).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

WB, IF, ChIP


## Dilutions

WB 1:500 - 1:2000

IF 1:50 - 1:200

ChIP 1:50 - 1:200

## Validations



Western blot - TEAD1 Polyclonal Antibody