

# TSHB Polyclonal Antibody



Catalog Number: 165038

## Product name

TSHB Polyclonal Antibody

## Specificity

Human, Mouse

## Antibody description

Polyclonal antibody to TSHB

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 21-138 of human TSHB (NP\_000540.2).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

WB, IHC, IF


## Dilutions

WB 1:500 - 1:2000

IHC 1:50 - 1:100


IF 1:50 - 1:100

## Validations



Immunohistochemistry - TSHB Polyclonal Antibody


Immunohistochemistry of paraffin-embedded human prostate using TSHB antibody at dilution of



1:100 (40x lens).

Immunohistochemistry - TSHB Polyclonal Antibody

Immunohistochemistry of paraffin-embedded mouse lung using TSHB antibody at dilution of



1:100 (40x lens).

Immunofluorescence - TSHB Polyclonal Antibody

Immunofluorescence analysis of U2OS cells using

# **TSHB Polyclonal Antibody**



Catalog Number: 165038

---

TSHB antibody . Blue: DAPI for nuclear staining.