

ZBED1 Polyclonal Antibody



Catalog Number: 165045

Product name

ZBED1 Polyclonal Antibody

Specificity

Human

Antibody description

Polyclonal antibody to ZBED1

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-280 of human ZBED1 (NP_001164607.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC, IF


Dilutions

WB 1:500 - 1:2000

IHC 1:50 - 1:200


IF 1:50 - 1:200

Validations



Western blot - ZBED1 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using ZBED1 Antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 60s.



Immunofluorescence - ZBED1 Polyclonal Antibody

Immunofluorescence analysis of U-2 OS cells using ZBED1 antibody at dilution of 1:100. Blue: DAPI for nuclear staining.