

SEPT7 Polyclonal Antibody

Catalog Number: 165070

Product name

SEPT7 Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to SEPT7

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 307-436 of human SEPT7 (NP_001011553.2).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC, IF


Dilutions

WB 1:500 - 1:2000

IHC 1:50 - 1:200


IF 1:50 - 1:100

Validations



Western blot - SEPT7 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using SEPT7 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 30s.




Immunohistochemistry - SEPT7 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human lung cancer using SEPT7 antibody at dilution of 1:100 (40x lens).


SEPT7 Polyclonal Antibody

Catalog Number: 165070



Immunohistochemistry - SEPT7 Polyclonal Antibody


Immunohistochemistry of paraffin-embedded mouse kidney using SEPT7 antibody at dilution of 1:100 (40x lens).



Immunofluorescence - SEPT7 Polyclonal Antibody
Immunofluorescence analysis of A-549 cells using SEPT7 antibody .

Immunohistochemistry - SEPT7 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded mouse Intestine using SEPT7 antibody at dilution



of 1:100 (40x lens).