

COL10A1 Polyclonal Antibody



Catalog Number: 165072

Product name

COL10A1 Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to COL10A1

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 521-680 of human COL10A1 (NP_000484.2).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

IHC

Dilutions


IHC 1:100 - 1:200

Validations



Immunohistochemistry - COL10A1 Polyclonal Antibody


Immunohistochemistry of paraffin-embedded rat lung using COL10A1 antibody at dilution of 1:100



(40x lens).

Immunohistochemistry - COL10A1 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human liver using COL10A1 antibody at dilution of



1:100 (40x lens).


Immunohistochemistry - COL10A1 Polyclonal Antibody

COL10A1 Polyclonal Antibody



Catalog Number: 165072

Immunohistochemistry of paraffin-embedded human stomach using COL10A1 antibody at dilution of 1:100 (40x lens).



Immunohistochemistry - COL10A1 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded mouse lung using COL10A1 antibody at dilution of 1:100 (40x lens).