

GABPB1 Polyclonal Antibody



Catalog Number: 165086

Product name

GABPB1 Polyclonal Antibody

Specificity

Human, Mouse

Antibody description

Polyclonal antibody to GABPB1

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 231-395 of human GABPB1 (NP_005245.2).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC, IF


Dilutions

WB 1:500 - 1:2000

IHC 1:50 - 1:200


IF 1:50 - 1:100

Validations




Western blot - GABPB1 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using GABPB1 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 90s.



Immunohistochemistry - GABPB1 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded mouse kidney using GABPB1 Antibody at dilution



of 1:100 (40x lens).


GABPB1 Polyclonal Antibody



Catalog Number: 165086


Immunohistochemistry - GABPB1 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human stomach using GABPB1 Antibody at dilution of 1:100 (40x lens).



Immunohistochemistry - GABPB1 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human gastric cancer using GABPB1 Antibody at dilution of 1:100 (40x lens).



Immunofluorescence - GABPB1 Polyclonal Antibody

Immunofluorescence analysis of MCF7 cells using GABPB1 antibody .