

HABP2 Polyclonal Antibody



Catalog Number: 165092

Product name

HABP2 Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to HABP2

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 24-313 of human HABP2 (NP_004123.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IF

Dilutions

WB 1:500 - 1:2000

IF 1:50 - 1:100


Validations



Western blot - HABP2 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using HABP2 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL

Basic Kit .Exposure time: 15s.



Immunofluorescence - HABP2 Polyclonal Antibody

Immunofluorescence analysis of MCF7 cells using HABP2 antibody . Blue: DAPI for nuclear staining.