

# HNRNPM Polyclonal Antibody



Catalog Number: 165105

## Product name

HNRNPM Polyclonal Antibody

## Specificity

Human, Mouse, Rat

## Antibody description

Polyclonal antibody to HNRNPM

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 491-730 of human HNRNPM (NP\_005959.2).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

WB, IHC, IF


## Dilutions

WB 1:500 - 1:2000

IHC 1:50 - 1:200


IF 1:50 - 1:200

## Validations




Western blot - HNRNPM Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using HNRNPM antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 60s.



Immunohistochemistry - HNRNPM Polyclonal Antibody

Immunohistochemistry of paraffin-embedded rat brain using HNRNPM antibody at dilution of 1:100




(40x lens).

Catalog Number: 165105

---


## Immunohistochemistry - HNRNPM Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human liver cancer using HNRNPM antibody at dilution of 1:100 (40x lens).




## Immunohistochemistry - HNRNPM Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human colon carcinoma using HNRNPM antibody at dilution of 1:100 (40x lens).




## Immunohistochemistry - HNRNPM Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human esophageal cancer using HNRNPM antibody at dilution of 1:100 (40x lens).



## Immunohistochemistry - HNRNPM Polyclonal Antibody


Immunohistochemistry of paraffin-embedded mouse kidney using HNRNPM antibody at dilution



of 1:100 (40x lens).

## Immunohistochemistry - HNRNPM Polyclonal Antibody

Immunohistochemistry of paraffin-embedded mouse heart using HNRNPM antibody at dilution of 1:100 (40x lens).



## Immunofluorescence - HNRNPM Polyclonal Antibody

# HNRNPM Polyclonal Antibody



Catalog Number: 165105

---

Immunofluorescence analysis of MCF7 cells using HNRNPM antibody .