

RPLP2 Polyclonal Antibody



Catalog Number: 165120

Product name

RPLP2 Polyclonal Antibody

Specificity

Human

Antibody description

Polyclonal antibody to RPLP2

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-115 of human RPLP2 (NP_000995.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

IF

Dilutions

IF 1:50 - 1:100

Validations



Immunofluorescence - RPLP2 Polyclonal Antibody

Immunofluorescence analysis of U2OS cells using RPLP2 antibody . Blue: DAPI for nuclear staining.