

NMU Polyclonal Antibody



Catalog Number: 165168

Product name

NMU Polyclonal Antibody

Specificity

Human, Mouse

Antibody description

Polyclonal antibody to NMU

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 35-174 of human NMU (NP_006672.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

IHC

Dilutions


IHC 1:50 - 1:100

Validations



Immunohistochemistry - NMU Polyclonal Antibody


Immunohistochemistry of paraffin-embedded human liver using NMU antibody at dilution of



1:100 (40x lens).

Immunohistochemistry - NMU Polyclonal Antibody


Immunohistochemistry of paraffin-embedded human esophagus using NMU antibody at dilution



of 1:100 (40x lens).

Immunohistochemistry - NMU Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human stomach using NMU antibody at dilution of



Immunohistochemistry - NMU Polyclonal Antibody

Immunohistochemistry of paraffin-embedded mouse lung using NMU antibody at dilution of 1:100 (40x lens).

1:100 (40x lens).