

# LMOD1 Polyclonal Antibody



Catalog Number: 165182

---

**Product name**

LMOD1 Polyclonal Antibody

**Specificity**

Human

**Antibody description**

Polyclonal antibody to LMOD1

**Preparation**

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 29-187 of human LMOD1 (NP\_036266.2).

**Formulation**

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Storage**

Store at -20°C. Avoid freeze / thaw cycles.

**Clonality**

Polyclonal

**Ig Type**

Rabbit IgG


**Applications**

WB, IHC

**Dilutions**

WB 1:500 - 1:2000

IHC 1:50 - 1:200

**Validations**

Immunohistochemistry - LMOD1 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human colon using LMOD1 antibody at dilution of 1:200 (40x lens).