

L3MBTL1 Polyclonal Antibody



Catalog Number: 165184

Product name

L3MBTL1 Polyclonal Antibody

Specificity

Human

Antibody description

Polyclonal antibody to L3MBTL1

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 533-772 of human L3MBTL1 (NP_056293.4).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

IF

Dilutions


IF 1:50 - 1:100

Validations



Immunofluorescence - L3MBTL1 Polyclonal Antibody

Immunofluorescence analysis of A549 cells using L3MBTL1 antibody .



Immunofluorescence - L3MBTL1 Polyclonal Antibody

Immunofluorescence analysis of GFP-RNF168 transgenic U2OS cells using L3MBTL1 antibody .Green:GFP-RNF168 fusion protein expression for DNA damage marker. Blue:DAPI for nuclear staining.RNF168 can be used to mark cells damaged by UV-A laser for they always gather around DNA damage region.