

CCAR2 Polyclonal Antibody



Catalog Number: 165201

Product name

CCAR2 Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to CCAR2

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 654-923 of human CCAR2 (NP_066997.3).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC, IF


Dilutions

WB 1:500 - 1:2000

IHC 1:50 - 1:200


IF 1:50 - 1:200

Validations




Western blot - CCAR2 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using CCAR2 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 30s.



Immunohistochemistry - CCAR2 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human kidney using CCAR2 antibody at dilution of




1:100 (40x lens).

CCAR2 Polyclonal Antibody


Catalog Number: 165201

Immunohistochemistry - CCAR2 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded rat brain using CCAR2 antibody at dilution of 1:100




(40x lens).



Immunohistochemistry - CCAR2 Polyclonal Antibody


Immunohistochemistry of paraffin-embedded mouse liver using CCAR2 antibody at dilution of



1:100 (40x lens).

Immunohistochemistry - CCAR2 Polyclonal Antibody


Immunohistochemistry of paraffin-embedded rat kidney using CCAR2 antibody at dilution of 1:100



(40x lens).

Immunohistochemistry - CCAR2 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded mouse brain using CCAR2 antibody at dilution of



1:100 (40x lens).

Immunohistochemistry - CCAR2 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human rectal cancer using CCAR2 antibody at dilution of 1:100 (40x lens).


Immunohistochemistry - CCAR2 Polyclonal Antibody

CCAR2 Polyclonal Antibody



Catalog Number: 165201

Immunohistochemistry of paraffin-embedded mouse kidney using CCAR2 antibody at dilution of



1:100 (40x lens).

Immunohistochemistry - CCAR2 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded mouse stomach using CCAR2 antibody at dilution of 1:100 (40x lens).