

# CSRNP1 Polyclonal Antibody



Catalog Number: 165205

## Product name

CSRNP1 Polyclonal Antibody

## Specificity

Human, Mouse, Rat

## Antibody description

Polyclonal antibody to CSRNP1

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-200 of human CSRNP1 (NP\_149016.2).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

WB, IHC, IF


## Dilutions

WB 1:500 - 1:2000

IHC 1:50 - 1:200


IF 1:50 - 1:100

## Validations



Immunohistochemistry - CSRNP1 Polyclonal Antibody


Immunohistochemistry of paraffin-embedded rat kidney using CSRNP1 antibody at dilution of 1:200



(40x lens).

Immunohistochemistry - CSRNP1 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded mouse brain using CSRNP1 antibody at dilution of 1:200 (40x lens).



Immunofluorescence - CSRNP1 Polyclonal Antibody

# CSRNP1 Polyclonal Antibody



Catalog Number: 165205

---

Immunofluorescence analysis of A549 cells using CSRNP1 antibody .