

NAP1L3 Polyclonal Antibody



Catalog Number: 165246

Product name

NAP1L3 Polyclonal Antibody

Specificity

Human

Antibody description

Polyclonal antibody to NAP1L3

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 307-506 of human NAP1L3 (NP_004529.2).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IF

Dilutions

WB 1:500 - 1:2000

IF 1:50 - 1:100

Validations



Immunofluorescence - NAP1L3 Polyclonal Antibody

Immunofluorescence analysis of A549 cells using NAP1L3 antibody . Blue: DAPI for nuclear staining.