

METTL7B Polyclonal Antibody

Catalog Number: 165248

Product name

METTL7B Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to METTL7B

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 55-244 of human METTL7B (NP_689850.2).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC, IF


Dilutions

WB 1:500 - 1:2000

IHC 1:50 - 1:200


IF 1:50 - 1:100

Validations




Western blot - METTL7B Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using METTL7B antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 90s.



Immunohistochemistry - METTL7B Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human stomach using METTL7B antibody at



dilution of 1:100 (40x lens).

METTL7B Polyclonal Antibody



Catalog Number: 165248

Immunofluorescence - METTL7B Polyclonal Antibody

Immunofluorescence analysis of A549 cells using METTL7B antibody . Blue: DAPI for nuclear staining.