

NAT8B Polyclonal Antibody



Catalog Number: 165250

Product name

NAT8B Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to NAT8B

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 75-165 of human NAT8B (NP_057431.2).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IHC

Dilutions

WB 1:500 - 1:2000


IHC 1:50 - 1:200

Validations




Western blot - NAT8B Polyclonal Antibody

Western blot analysis of extracts of rat brain, using NAT8B antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 90s.



Immunohistochemistry - NAT8B Polyclonal Antibody

Immunohistochemistry of paraffin-embedded rat kidney using NAT8B antibody at dilution of 1:200



(40x lens).


NAT8B Polyclonal Antibody



Catalog Number: 165250

Immunohistochemistry - NAT8B Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human stomach using NAT8B antibody at dilution



of 1:200 (40x lens).

Immunohistochemistry - NAT8B Polyclonal Antibody

Immunohistochemistry of paraffin-embedded mouse heart using NAT8B antibody at dilution of 1:200 (40x lens).