

# SNRPD1 Polyclonal Antibody



Catalog Number: 165254

---

## Product name

SNRPD1 Polyclonal Antibody

## Specificity

Human

## Antibody description

Polyclonal antibody to SNRPD1

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-119 of human SNRPD1 (NP\_008869.1).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

IF

## Dilutions

IF 1:50 - 1:100

## Validations



Immunofluorescence - SNRPD1 Polyclonal Antibody

Immunofluorescence analysis of MCF-7 cells using SNRPD1 antibody .