

TAF5 Polyclonal Antibody



Catalog Number: 165258

Product name

TAF5 Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to TAF5

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 656-800 of human TAF5 (NP_008882.2).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IF

Dilutions

WB 1:500 - 1:2000


IF 1:50 - 1:100

Validations



Western blot - TAF5 Polyclonal Antibody

Western blot analysis of extracts of mouse brain, using TAF5 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time:



Immunofluorescence - TAF5 Polyclonal Antibody

Immunofluorescence analysis of A-549 cells using TAF5 antibody. Blue: DAPI for nuclear staining.