

LMO4 Polyclonal Antibody



Catalog Number: 165283

Product name

LMO4 Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to LMO4

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-165 of human LMO4 (NP_006760.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


IHC, IF

Dilutions

IHC 1:50 - 1:200


IF 1:50 - 1:100

Validations



Immunohistochemistry - LMO4 Polyclonal Antibody


Immunohistochemistry of paraffin-embedded rat brain using LMO4 antibody at dilution of 1:100



(40x lens).

Immunohistochemistry - LMO4 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human placenta using LMO4 antibody at dilution



of 1:100 (40x lens).


Immunohistochemistry - LMO4 Polyclonal Antibody

LMO4 Polyclonal Antibody



Catalog Number: 165283

Immunohistochemistry of paraffin-embedded mouse heart using LMO4 antibody at dilution of



1:100 (40x lens).

Immunofluorescence - LMO4 Polyclonal Antibody

Immunofluorescence analysis of MCF-7 cells using LMO4 antibody . Blue: DAPI for nuclear staining.