

ATOH7 Polyclonal Antibody



Catalog Number: 165287

Product name

ATOH7 Polyclonal Antibody

Specificity

Human

Antibody description

Polyclonal antibody to ATOH7

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 47-152 of human ATOH7 (NP_660161.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IF

Dilutions

WB 1:500 - 1:2000

IF 1:50 - 1:100

Validations



Immunofluorescence - ATOH7 Polyclonal Antibody

Immunofluorescence analysis of A-549 cells using ATOH7 antibody .