

GCM1 Polyclonal Antibody



Catalog Number: 165289

Product name

GCM1 Polyclonal Antibody

Specificity

Human

Antibody description

Polyclonal antibody to GCM1

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 207-436 of human GCM1 (NP_003634.2).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IF

Dilutions

WB 1:500 - 1:2000

IF 1:50 - 1:100

Validations



Immunofluorescence - GCM1 Polyclonal Antibody

Immunofluorescence analysis of MCF-7 cells using GCM1 antibody . Blue: DAPI for nuclear staining.