

RNF7 Polyclonal Antibody



Catalog Number: 165298

Product name

RNF7 Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to RNF7

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-113 of human RNF7 (NP_055060.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IF

Dilutions

WB 1:500 - 1:2000

IF 1:20 - 1:50

Validations



Immunofluorescence - RNF7 Polyclonal Antibody

Immunofluorescence analysis of MCF-7 cells using RNF7 antibody . Blue: DAPI for nuclear staining.