

TFEB Polyclonal Antibody

Catalog Number: 165304

Product name

TFEB Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to TFEB

Preparation

Antigen: Recombinant protein of human TFEB.

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC, IF


Dilutions

WB 1:500 - 1:2000

IHC 1:50 - 1:200


IF 1:50 - 1:200

Validations




Western blot - TFEB Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using TFEB antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 30s.



Immunofluorescence - TFEB Polyclonal Antibody

Immunofluorescence analysis of HeLa cells using TFEB antibody at dilution of 1:100. Blue: DAPI for



nuclear staining.

TFEB Polyclonal Antibody



Catalog Number: 165304

Immunofluorescence - TFEB Polyclonal Antibody

using TFEB antibody at dilution of 1:100. Blue:
DAPI for nuclear staining.

Immunofluorescence analysis of NIH/3T3 cells