

PTGS1 Polyclonal Antibody



Catalog Number: 165318

Product name

PTGS1 Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to PTGS1

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-180 of human PTGS1 (NP_001258094.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC, IF


Dilutions

WB 1:500 - 1:2000

IHC 1:100 - 1:200


IF 1:50 - 1:200

Validations



Western blot - PTGS1 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using PTGS1 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure



time: 15s.

Immunohistochemistry - PTGS1 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human liver cancer using PTGS1 antibody at dilution of 1:100 (40x lens).