

PEX3 Polyclonal Antibody



Catalog Number: 165324

Product name

PEX3 Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to PEX3

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 144-373 of human PEX3 (NP_003621.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IHC

Dilutions

WB 1:500 - 1:2000


IHC 1:50 - 1:200

Validations



Western blot - PEX3 Polyclonal Antibody


Western blot analysis of extracts of various cell lines, using PEX3 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure



time: 90s.

Immunohistochemistry - PEX3 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human liver injury using PEX3 antibody at dilution




of 1:100 (40x lens).

PEX3 Polyclonal Antibody

Catalog Number: 165324

Immunohistochemistry - PEX3 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded rat kidney using PEX3 antibody at dilution of 1:100



(40x lens).

Immunohistochemistry - PEX3 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded mouse kidney using PEX3 antibody at dilution of 1:100 (40x lens).