

SH3GL1 Polyclonal Antibody



Catalog Number: 165371

Product name

SH3GL1 Polyclonal Antibody

Specificity

Human, Mouse, Rat, Monkey

Antibody description

Polyclonal antibody to SH3GL1

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 219-368 of human SH3GL1 (NP_003016.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IHC

Dilutions

WB 1:500 - 1:2000

IHC 1:100 - 1:200

Validations



Western blot - SH3GL1 Polyclonal Antibody

Western blot analysis of extracts of COS-7 cells, using SH3GL1 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 60s.