

# TTBK2 Polyclonal Antibody



Catalog Number: 165454

## Product name

TTBK2 Polyclonal Antibody

## Specificity

Human, Rat

## Antibody description

Polyclonal antibody to TTBK2

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1185-1244 of human TTBK2 (NP\_775771.3).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


WB, IHC

## Dilutions

WB 1:500 - 1:2000


IHC 1:50 - 1:200

## Validations



Western blot - TTBK2 Polyclonal Antibody

Western blot analysis of extracts of rat brain, using TTBK2 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 90s.



Immunohistochemistry - TTBK2 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human stomach using TTBK2 antibody at dilution of 1:200 (40x lens).