

ACE Polyclonal Antibody



Catalog Number: 165479

Product name

ACE Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to ACE

Preparation

Antigen: A synthetic peptide of human ACE

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IHC

Dilutions

WB 1:500 - 1:2000


IHC 1:50 - 1:100

Validations




Western blot - ACE Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using ACE antibody at 1:500 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 90s.



Immunohistochemistry - ACE Polyclonal Antibody

Immunohistochemistry of paraffin-embedded mouse brain using ACE antibody (40x lens).



ACE Polyclonal Antibody



Catalog Number: 165479

Immunohistochemistry - ACE Polyclonal Antibody

Immunohistochemistry - ACE Polyclonal Antibody

Immunohistochemistry of paraffin-embedded rat brain using ACE antibody (40x lens).

Immunohistochemistry of paraffin-embedded human esophagus using ACE antibody (40x lens).

