

SOCS1 Polyclonal Antibody



Catalog Number: 165519

General Information

Product name:

SOCS1 Polyclonal Antibody

Immunogen:

Applications:

WB, IHC

Specificity:

Human, Mouse, Rat

Antibody description:

Polyclonal antibody to SOCS1

Properties

Preparation:

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-211 of human SOCS1 (NP_003736.1).

Formulation:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage:

Store at -20°C. Avoid freeze / thaw cycles.

Clonality:

Polyclonal


Ig Type:Rabbit IgG

Applications

WB 1:500 - 1:2000

IHC 1:50 - 1:200

Validations



Western blot - SOCS1 Polyclonal Antibody

Western blot analysis of extracts of mouse lung, using SOCS1 antibody at 1:3000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 60s.